

## STEP 6

**Activity Wristband navigation.** To navigate the Activity Wristband use the touch key located at the bottom of the display. Short pressing the touch key moves from one screen to another. Long pressing the touch key enters and exits submenus.

For example (accessing the Activity screen):

- To move from Home screen to Activity screen, short press the touch key.
- On Activity long press the touch key which enters the activity submenus starting with the Steps.
- Short press the touch key and move from Steps to Distance to Calories to BPM to Back to Steps and so on.
- To exit the Activity submenu, navigate to Back and long press the touch key. You are now back on the main Activity screen.

There are 9 main screens and several submenus. The 9 main screens are: Home, Activity, Sleep, Sport, Notifications, Weather, Alarm, Tools, Settings.

## Warranty Registration

Visit us online to register your product and receive special offers and news from Zewa.

For a detailed description of your warranty, please refer to the instruction manual.

[www.zewa.com](http://www.zewa.com)



Zewa Inc.  
12960 Commerce Lakes Dr. # 29  
Fort Myers, FL 33913  
Toll Free: 888-993-3592

[www.zewa.com](http://www.zewa.com)



## Getting Started

Warranty Registration



Activity Wristband  
(Ver 3)

## STEP 1



### Charge Activity Wristband before pairing it to the App.

Remove wristband (with holes) from the Activity Wristband

NOTE: Removing the band will get easier over time.

## STEP 2

Use your phone charger to charge the Activity Wristband (charger not included.)

To use your phone charger, remove cable from the wall charging unit and insert the Activity Wristband.

Charge for about 5 hours.



Note 1: The white part inside the wall charger needs to be at the bottom.

Insert Activity Wristband, display side up.

Note 2: If the Activity Wristband does not charge after 15sec. remove

Activity Wristband, flip charger around and make sure the white part inside charger is at the bottom.

Insert Activity Wristband, display side up.

## STEP 3

Download and install the Zewa-Fit App at [www.zewa.com/zewa-fit](http://www.zewa.com/zewa-fit)

IMPORTANT: During the install (or when first opening the app), the app will ask several installation questions, accept (say "yes") to all of them.

## STEP 4

Create a new user account in the App.

- After opening the App for the first time, login using an existing login option or create a new user account by selecting "Other Login Options".
- When selecting "Other Login Options", select "Create Account".
- Enter a valid E-Mail address and press "Get Code". The system will E-Mail you a verification code. Enter that code into the app and press "Login".
- Continue with entering your information.

## STEP 5

Pair Activity Wristband to the App.

- In the app, tap on "Device" on the bottom of your screen.
- Tap on "Add device".
- Select device you would like to add.
- Confirm pairing by pressing the touch key on the Activity Tracker.