



User Manual

Bathroom Scale

Model # 21301



Manufactured for Zewa, Inc.
12960 Commerce Lakes Drive # 29
Fort Myers, FL 33913
USA
www.zewa.com
Toll Free Customer Service:
1-888-993-3592
warranty@zewa.com







- Thank you very much for selecting the Zewa Bathroom Scale 21301.
- Please do read the user manual carefully and thoroughly so as to ensure the safe usage of this product, and keep the manual well for further reference in case you have problems.

Table of Contents

Safety Information	
Safety and Usage Information	2
Your Scale and Its Environment	3
Efficient Use of Your Scale	3
Overview	
General Instructions	4
Device Components	4
Function Indication	4
List	5
LCD Display	5
Initial Start-Up	
Insert the Batteries	6
Start Measuring	
Select Measurement Unit	7
Initialising Your Scale.....	8
Start Measuring	8
Troubleshooting	
Error Prompt	9
Specifications	10
Maintenance	11
Warranty	11
Appendix	
FCC Regulations.....	12
EMC Guidance	12

♥ **Safety and Usage Information**

The warning signs and symbols are essential to ensure your correct and safe use of this product and protect you and others from injury. Please kindly find the meanings of the warning signs and symbols, which you may encounter in the label and user manual, as follows:

	Symbol for "THE OPERATION GUIDE MUST BE READ"	SN	Symbol for "SERIAL NUMBER"
	Symbol for "DIRECT CURRENT"		Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "COMPLIES WITH EU REQUIREMENTS"		

♥ **Your Scale and Its Environment**

To ensure your safety as well as the service life of your scale, please avoid using the scale under the following circumstances:

- Slippery floor such as tile floor
- Jumping onto the platform immediately after bath or with wet feet
- Near a cell phone or microwave oven

Avoid storage in the following locations:

- Where there is water
- Where the device may be exposed to extreme temperatures, humidity, moisture, direct sunlight, dust, or salt air
- Where there is risk of shock or drop
- Where you store chemicals or full of corrosive gases
- Where in reach of the infants or children

♥ **Efficient Use of Your Scale**


To ensure the accuracy of measurement, please follow below instructions when you start measurement.

- Place the scale on a flat, hard surface. Soft surface such as carpet will affect the performance of the scale.
- Start measurement at least two hours after Getting up or Dinning.
- Avoid measurement immediately after strenuous exercise, sauna or bath, drinking, and dinning.
- Always start measurement in the same time slot and on the same scale located on the same flat, hard surface.
- In order to ensure the accuracy of weighing, when you moved the scale or haven't used the scale for a long time, please turn on the scale first, wait until it turns off, then you can stand on the platform again and measure your weight.
- When you don't use the scale, please place it stable. Do not put anything on it.
- Any information provided by this scale is in no way meant to treat, cure or prevent any disease or illness from happening. If in doubt, contact your physician.
- It is intended for use in the domestic setting only.

♥ **General Instructions**

Zewa Bathroom Scale 21301 applies SENSE ON technology. It's not necessary to turn on or turn off the scale, just stand on the scale for measurement. It is delicate and trendy, convenience for using.

♥ **Device Components**




♥ **Function Indication**

◇ UNIT Button: Press "UNIT Button" to select desired unit in weighing mode.
or
Press "UNIT Button" to turn on the scale.

♥ **List**

1. Bathroom Scale 21301
2. One CR2032 Battery (3V)
3. User Manual

♥ **LCD Display**



kg	Kilogram
lb	Pound
st	Stone
lb	

♥ Insert the Battery

- Open the battery door in the back of the scale.
 - Insert the battery (1 x CR2032) into the battery compartment according to the polarity indications marked inside the compartment.
- * All the digits will be shown on the LCD.



- Close the battery door and wait until the digits "0.0lb" are shown on the LCD.




CAUTION

- When the symbol $\downarrow 0$ appears, the device will power off in four seconds. Then you shall replace with a new one battery.
- Worn batteries are hazardous waste. Do NOT dispose of them together with the household garbage. Please refer to the local ordinances and recycling instructions regarding disposal of the worn batteries.

♥ Select Measurement Unit

With battery correctly installed, press "UNIT" button in the back of the scale to select measurement unit. The default measurement unit is "pound". You may press "UNIT" button to choose among pound, kilogram and stone.



CAUTION

Please put the scale on the flat hard ground when finishing the unit selection. Wait for a few seconds, then you may start measuring according to the measuring procedure below. If there is not any operation in weighing mode, the scale will turn off automatically after about 10s.

♥ Initialising Your Scale

1. Press the platform centre and remove your foot
or
Press UNIT button to turn on the scale.
2. "0.0" will be displayed.



3. The scale will switch off and is now ready for use.
This initialisation process must be repeated if the scale is moved.
At all other times step straight on the scale.

♥ Start Measuring



1. Place the scale flatly on a hard surface.
2. Step onto the platform and stand still until the numeral on the LCD stops blinking, then the weight is locked.



3. The measuring result is displayed for about 10 seconds, and then it will turn off.

Note: If there is no operation in 10 seconds when the scale is on, it will turn off automatically.

♥ Error Prompt

Error	Description	Solution
	Overload. The device will power off.	Stop using this scale for measurement.
	Low Battery. The device will power off in four seconds.	Replace with 1xCR2032 new battery. Please purchase the authorized batteries for replacement.

♥ Specifications

Product Name	Zewa Bathroom Scale 21301
Dimension	Scale: 300x300x19mm (Approximately)
Net Weight	Approximately 1.52kg (Excluding the cell)
Display	Digital LCD
Measurement Unit	Pound / Kilogram/ Stone
Measurement Range	330.6 lb / 150kg / 23 st: 8 lb
Division	0.1kg / 0.2 lb/1 lb (≥20st)
Accuracy	50kg: ±0.3kg; 100kg: ±0.4kg; 150kg: ±0.5kg;
Working Environment	Temperature: 5 C to 40 C Relative Humidity: ≤90% RH Atmospheric pressure: 86kPa to 106kPa
Storage Environment	Temperature: -20 C to 60 C Relative Humidity: ≤90% RH Atmospheric pressure: 50kPa to 106kPa
Power Source	3V (One CR2032 Battery)
Auto-ON	SENSE ON technology
Auto-OFF	The scale will turn off after about 10s if there is no operation.
Accessories	1. One CR2032 Battery 2. User Manual
Mode of Operation	Continuous Operation

About the Accuracy of This Product

- This product passes strict inspection before delivery and therefore its accuracy is guaranteed by the manufacturer. Please refer to the above table for the descriptions on accuracy.
- This product is specially designed for body scale as well as weight measurement. It should NOT be used by anyone during the process of transaction for verification of goods' weight.

♥ Maintenance

When carrying out usual maintenance, please ensure practice of the following Do's and Don'ts:

- DO use a dry soft cloth to wipe the dust.
- DO use a wet soft cloth, dipped into water and wrung out, to wipe the dirt. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use propellant, abrasive or other chemicals to wipe the dirt in avoidance of discolor or malfunction.
- DON'T disassemble this device. If you have any problems, please contact Zewa.
(Please refer to Warranty for contact information)

♥ Warranty

- Zewa warrants its products free of defects in materials and workmanship in normal use for a period of TWO years from the date of retail purchase.
- This warranty does NOT cover damages caused by misuse or abuse, including but not limited to:
 - Failure caused by unauthorized repairs or modifications;
 - Damage caused by shock or drop during transportation;
 - Failure caused by improper operation inconsistent with the instructions stated in this user manual;
 - Malfunction or damage from failure to provide the recommended maintenance;
 - Damage caused by improper use of power supply.
- Should this device require maintenance (or replacement at our option) under warranty, please deliver the original package to Zewa Inc. prepaid. Please return the store receipt (with the retail purchase date) and a note with reasons to return on it as well.

Zewa, Inc.
12960 Commerce Lakes Drive # 29, Fort Myers, FL 33913, USA
www.zewa.com
Toll Free Customer Service: 1-888-993-3592
warranty@zewa.com

♥ FCC Regulations

FCC User Guide Information

Radio Frequency Interface Requirements - FCC

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna;
- Increase the separation between the equipment and receiver;
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected;
- Consult the dealer or an experienced radio / TV technician for help.

Radio Transmitters (Part 15)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

♥ EMC Guidance

1. MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
2. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d = 3,3$ m away from the equipment. (Note. As indicated in Table 6 of IEC 60601-1-2:2007 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields $d = 3,3$ m at an IMMUNITY LEVEL of 3 V/m)