



Patient Onboarding Instructions

For any issues or additional questions, please contact Zewa at 1-888-993-3592 / 239-337-9975 or support@zewa.com.



Online Patient Set Up

www.myhealthconnected.net / Log in with your credentials

To add a new patient, click on “Dashboard” - “Patient Intake Form”.

The screenshot shows the MyHealth Connected dashboard for a user named Zewa - Customers. The dashboard includes a navigation bar with links for DASHBOARD, PATIENTS, ADMIN, ABOUT, and LOG OUT. A red box highlights the PATIENT INTAKE FORM button in the MESSAGE CENTER (WIP) section. The main content area displays a welcome message, user profile information, and various statistics for the user's population, including compliance rates and data posting status.

Enter patient information, then click on “Add Patient”.

The screenshot shows the Patient Intake Form in the MyHealth Connected system. The form is titled "Patient Intake Form" and includes fields for patient information: Email, First Name, Last Name, Middle Initial, Measurement System (US/Imperial), Gender & Weight, Height (Feet and Inches), DOB (Month, Day, Year), and Time Zone. A red box highlights the "Add Patient" button at the bottom of the form.

The next screen is used to add devices to the patient. Select “Blood Pressure Monitor” from the drop down menu and then add the serial number of the device. Each device has a sticker that shows a serial number and mac address, use the serial number on that sticker (sample sticker shown below). When finished, click on “Add New Device”.



My Health Connected

DASHBOARD PATIENTS ADMIN ABOUT LOG OUT

Zeindler

Edit Devices

Action	Device Type	Identifier	Entry Date
There are not any devices assigned...			

Add a new device

Device Type Blood Pressure Monitor (Zewa)

User# User 1

Serial#

Click 2 times → Add New Device Cancel Close

The newly added device is listed in the device list. To add more devices, repeat this step, otherwise click on close.

Now we are back on the Patient Intake Form. Click on “Go to Patient Dashboard”, to finish the patient and device set up.

NOTE: The patient will receive a system generated email to verify this account. They only have to follow the E-mail instructions IF they want to access their own information online. Otherwise they can ignore this E-Mail and no further action is needed.

This last screen is VERY IMPORTANT. For app set up, please note the E-Mail address and Password listed below. Leave this page open and proceed with the app set up on the patients mobile device.

Android, Test

Android, Test

Android, Test

Add/Edit Devices

On Hold

Gender: Male, Date of Birth: 1/1/1970

Patient Local Time 11:52 AM -05:00

Last portal login: No previous login recorded...

Last Data Import: No device contact yet...

Email: demo19@zewa.com

Temporary Password: L86684xx

Message Center: 0 New Messages

Blood Pressure

Weight

Activity

Last Data: N/A

Last Data: N/A

Last Data: N/A

Systolic

Diastolic

Pulse

P-Press.

Current

Previous Var.

Base Var.

Steps

kCal

Last 7 Day's Readings

Last 7 Day's Readings

Last 7 Day's Readings

Average

Maximum

Minimum

Average

Maximum

Minimum

Sum Var.

Average


Maximum

Minimum

App Set Up for Android Devices

On the mobile device where the app has to be installed, confirm that Bluetooth is turned ON and then open the Google “Play Store”. Search for “myhealthconnected” or “zewa” and select the “**MyHealthConnected Zewa Inc**” app (the app symbol is a red heart with a stethoscope wrapped around it). Install and open the app. Then follow these instructions:

- Click on “LOG IN”.
- Enter login and password from online portal (below is just a sample, do NOT copy this). If you left your screen open as recommended, copy information.



- A window will pop up saying that “No Measurement Devices Defined”. Click on “Yes” to add a new device.
- Checkmark “BP Monitor” and click on “Done”. Follow pairing instructions.
- Once pairing is finished, click on “OK” and then “Done”.


Your set up is now complete, please train patient on using the Blood Pressure Monitor correctly.

NOTE: The app will collect data in the background as long as Bluetooth is left ON (on the mobile device). For a manual sync, open the app and tap on the “START SYNC” button.

App Set Up for iOS (Apple) Devices

On the mobile device where the app has to be installed, confirm that Bluetooth is turned ON and then open the Apple “App Store”. Search for “myhealthconnected” or “zewa” and select the “**zewa - myhealthconnected**” app (the app symbol is a red heart with a stethoscope wrapped around it. Install and open the app. Then follow these instructions.

- Click on “JOIN MyHealthConnected”.
- Follow the registration process. For E-Mail and password use the same credentials as in the online portal (below is just a sample, do NOT copy this). If you left your screen open as recommended, copy information.



- Click on “Finish Registration”.
- Once the app is open, tap on the “wheel symbol” in the top right corner.
- In the new menu select “measurement device”.
- Tap on the “Connect” button next to “Blood Pressure Monitor” and follow the instructions.
- Once pairing is finished, the device dashboard shows the word “Paired” next to “Blood Pressure Monitor”.
- Click the back arrow in the top left corner 2 times to get back to the main screen.

Your set up is now complete, please train patient on using the Blood Pressure Monitor correctly.

NOTE: The app will only collect data from your Bluetooth device if the app is open and in the foreground. Please open your app periodically, with Bluetooth turned ON and your medical device close by. This will connect to your medical device(s) and download the latest measurements.