



Cellular Blood Pressure Monitoring

Cellular BP Monitor | Model No. UAM-950

- Extra Large Display with Backlight
- 60 Reading Memory
- Irregular Heartbeat Indicator
- Date and Time Function
- 4G Cellular (AT&T Network)
- Wide Range Cuff: 8.6"-16.5"
- Batteries Included
- Compatible with:
Zewa-Connect Device Management Portal



*Designed exclusively for
Remote Patient Monitoring*



60 Reading Memory:

Monitors and tracks 60 readings for one user.



Cellular Enabled

4G Cellular (AT&T Network)



Irregular Heartbeat Indicator:

If irregular heartbeat is detected, an indicator icon will appear alerting you so you can consult your medical professional.



Zewa Connect

Compatible with Zewa's Proprietary Device Management Portal. An exclusive Zewa solution for Connected Devices.



Date and Time Function:

Automatically stores the date and time when a blood pressure reading is taken.