



## Infrared Thermometer | Rotating Sensor



Patient Measurement



Off Position



Patient Measurement

### Thermometer | Model No. 11110<sup>v2</sup>

- Forehead Non-Touch Technology
- Innovative 4-Way Rotating Temperature Sensor
- Distance Detection for Reliable, Accurate Readings
- LCD Display with 3-Color Temperature Guidance
- Fast 1-Second Measurement
- 30-Reading Memory Recall
- Two Mode Measurements, Person or Objects
- Measures in Celsius or Fahrenheit
- Power Supply: 2 x AAA Batteries.
- 2 Year Warranty



4-Way Rotating Temperature Sensor with Distance Detection Feature ensures Precise Readings

### Packaging and Shipping Specifications

	Dimensions (L x D x H)	Weight	Cube	UPC/GTIN
Selling Unit (Box)	1.8 x 4.3 x 9	0.4 lbs.	0.04	6-82891-11110-4
Selling Unit (10 Units)	9 x 9.5 x 9.5	4.4 lbs.	.47	1-0682891-11110-1
Master Case (40 Units)	11 x 18.8 x 20	19 lbs.	2.39	2-0682891-11110-8

Zewa Inc. | 12960 Commerce Lakes Dr. # 29 | Fort Myers, FL 33913

1-888-993-3592 | www.zewa.com