



Quality products for life



## Ear - Forehead Thermometer | Bluetooth



### Thermometer | Model No. 11115

- Ear or Forehead **Non-Touch** technology
- Professional accuracy measurement for all ages.
- Dual mode to measure via forehead and ear.
- Fast measurement with fever alarm.
- Measures in Celsius and Fahrenheit.
- Bluetooth connectivity.
- Power Supply: 2 x AAA Batteries.
- 1 Year Warranty
- UPC: 6-82891-11115-9



Infrared Thermometer measures the body temperature based on the infrared energy emitted from the ear or forehead. Users can quickly get measurement results after positioning the temperature probe at ear or forehead.

## Packaging and Shipping Specifications

	Dimensions (L x D x H)	Weight	Cube	UPC/GTIN
Selling Unit (Box)	1.8 x 1.9 x 6.1	0.3 lbs	0.1	6-82891-11115-9
Selling Unit (10 Units)	12 x 9 x 4	3.1 lbs	0.3	1-0682891-11115-6
Master Case (50 Units)	11 x 14.2 x 20.5	29.8 lbs	1.9	2-0682891-11115-3

Zewa Inc. | 12960 Commerce Lakes Dr. # 29 | Fort Myers, FL 33913

1-888-993-3592 | www.zewa.com