

# Healthy Living



## Wireless BMI Smart Scale



### Model | 21300

- 12 Body Composition Metrics
- Large, easy to read digital display
- Safe, Low Profile Design (3/4")
- Measures in lbs. and Kg
  - \*Up to 397 lbs. / 180 Kg
- Strong ITO Safety Glass
  - \* (Indium tin oxide coated glass)
- Bluetooth 5.0
- UPC # 6-82891-21300-6



**Scale measures in lbs and Kg** for up to 397 lbs / 180 Kg.



#### ITO Safety Glass

Sensors to measure body composition are integrated in the safety glass.



**Measure up to 12 body composition metrics**  
Includes BMI, BMR, body fat rate, muscle rate, water rate, body fat mass, muscle mass, bone mass, fat-free body weight and heart rate (calculated in app).



#### Wireless Data Transfer


Connects to Android and iOS mobile devices.



#### Take Control

Add additional Zewa Bluetooth devices to monitor more vital signs using the same app.

## Bluetooth Connectivity



### Measures up to 12 Body Composition Metrics

Includes BMI, BMR, body fat rate, muscle rate, bone rate, water rate, body fat mass, muscle mass, fat free body weight, and heart rate.

## Technical Specifications

<b>Product Dimensions</b>	<b>310mm x 310mm x 25mm (12.2" x 12.2" x 1")</b>
<b>Product Weight</b>	1588g / 3.5 lbs (w/o batteries)
<b>Display Size</b>	70mm / 2.8"
<b>Power Source</b>	4 x AAA Batteries

<b>Display</b>	<b>Digital display with Blue backlight</b>
<b>Measuring range</b>	5 - 180 Kg / 11 - 397 lbs
<b>Turn ON method</b>	Step ON sensor
<b>Accuracy</b>	+/- 3 mmHg

## Packaging and Shipping Specifications

	<b>Dimensions (L x D x H)</b>	<b>Weight</b>	<b>Cube</b>	<b>UPC/GTIN</b>
<b>Selling Unit (Box)</b>	13.5 x 13.5 x 2	4.5 lbs	0.21	6-82891-21300-6
<b>Inner Case (4 Units)</b>	14 x 14 x 8	18.7 lbs	0.91	1-0682891-21300-3
<b>Master Case (8 Units)</b>	15 x 14.5 x 17	39.4 lbs	2.14	2-0682891-21300-0