



Compressor Nebulizer System | Aerosol Delivery

Nebulizer | Model No. NEB-36555



- Aerosol Delivery for Patients of all Ages
- Intended for the Treatment of Asthma, COPD, Allergies and other Respiratory Disorders
- High-Performance, Insurance Compliant System
- Intuitive Design is Easy to Use and Maintain
- Durable Quiet Piston Driven System
- Includes Complete Patient Administration Set, 6 Air Filters
- Compact, Lightweight Design with Integrated Handle
- Consistently Delivers Medication ≤ 5 microns MMAD
- Provides Fast, Efficient Treatment Times
- 5 Year Limited Warranty



Advanced System exceeds market standards, offering superior results and performance

CONTENTS | Model NEB-36555 Single Patient Use Compressor Nebulizer System, Includes Compressor Nebulizer Unit, Medication Cup, Nebulizer Kit, Inhalation Tube, Corrugated Pipe, T-Type Connector, Inhalation Mouthpiece (1 ea.), Filter (6 pcs.) - [8 units per Case pack].

HCPCS Compliant E0570

Zewa Inc. | 12960 Commerce Lakes Dr. # 29 | Fort Myers, FL 33913

1-888-993-3592 | www.zewa.com

ZEWA Nebulizer | Model No. NEB-36555

ORDERING INFORMATION

MFR #	DESCRIPTION	CONTENTS	QUANTITY
NEB-36555	Compressor Nebulizer System	Compressor, Patient Neb Kit, *(Medication-Cup, Kink Resistant Air Tubing, Inhalation Corrugated Pipe, T-Type Connector, Mouthpiece), 6 Air-filters	1 Each (8 per case)
36556FLT	Nebulizer Filters	Air filters for NEB-36555	10/Pack (20 pks per case)
36557NKT	Patient Administration Kit *(Neb-Kit)	Medication Cup, Air Tubing, Inhalation Corrugated Pipe, Mouthpiece, T-Type Connector,	1 Each (20 per case)



SATISFACTION GUARANTEED

For more information regarding Zewa Brand products please visit us www.zewa.com

Technical Specifications

Product Dimensions	6.9 x 5 x 4
Product Weight	2.7 lbs.
Medication Capacity	5ml
Particle size	≤ 5 μm MMAD
Warranty	5 Year Limited Warranty

Max Compressor Flow	> 7 LPM
Max Pressure - Low Pressure	30 PSI
Operating Pressure Range	8-16 psi
Operating Noise Level	≤ 55 db

Packaging and Shipping Specifications

Product Packaging	Dimensions (L x D x H)	Weight	UPC/GTIN
Selling Unit (1 Box)	8.5 x 8.5 x 6	3 lbs.	6-82891-36555-2
Master Case (8 Units)	17.7 x 11.4 x 13.2	24 lbs.	2-0682891-36555-6

Zewa Inc. | 12960 Commerce Lakes Dr. # 29 | Fort Myers, FL 33913

1-888-993-3592 | www.zewa.com